




Protective Padded Helmets, Soft Shell Z9820

ITEM #	Z9820-[size]
DESCRIPTION	Protective Padded Helmets, Soft Shell
SIZE	L-XL (see next page for helmet sizing information)
SPECIFICATIONS	<ul style="list-style-type: none">• Helps reduce head Injuries• Feature ventilation holes and adjustable chin straps• Excellent head protection for transport
CASE PACK	EA
OTHER	This product is not a substitute for careful and regular monitoring of the individuals who utilize it, and this product must be used in strict conjunction with all laws, rules, regulations, policies and procedures that apply to the supervision, monitoring, and observation of persons who are or who are believed to be at risk of committing self-harm or suicide.

Helmet Sizing Information



***IMPORTANT: All measurements are in inches.**

(A) Circumference Range (all the way around the head)		
Size	Soft Shell	Hard Shell
Infant	16 1/4 - 17	Not available.
XX-Small	17 1/2 - 18 1/4	18 1/2 - 19 1/4
X-Small	19 - 19 3/4	19 1/4 - 20
Small	20 3/4 - 21 1/2	20 1/2 - 21 1/4
Medium	21 1/2 - 22 1/4	21 1/2 - 22 1/4
Large	23 - 23 3/4	22 1/2 - 23 1/4
X-Large	24 - 24 3/4	24 1/2 - 25 1/4
XX-Large	25 1/2 - 26 1/4	26 - 26 3/4

(B) Occipital		
Size	Soft Shell	Hard Shell
Infant	11 1/2	Not available.
XX-Small	12 3/4	12 3/4
X-Small	13 3/4	13 1/4
Small	14 1/2	15
Medium	15 1/4	16 1/2
Large	16 1/2	16 3/4
X-Large	17 1/2	17 1/2
XX-Large	18 1/2	18 3/4

(C) Ear To Ear		
Size	Soft Shell	Hard Shell
Infant	8 1/4	Not available.
XX-Small	9 1/2	9 1/4
X-Small	9 1/2	10
Small	10 1/2	11 3/4
Medium	11	12 1/2
Large	12	12 3/4
X-Large	13	13 1/4
XX-Large	13 3/4	15

Since proper fit is important for the comfort of the person wearing headgear, measurements should be determined by professional personnel.

****PLEASE MEASURE IN INCHES ONLY****

Measure the head circumference (A) at eyebrow level all the way around the head.

Occipital measurement (B) is taken from eyebrow level to the back of the head at the point that the helmet is to end.

Over the top of the Head (C) from the tip of one ear to the tip of the other ear.

Helmets should fit snug- special sizing pads are available for purchase. *See helmet accessories and replacement items for headgear.

If for some reason the measurements do not fall within the above guidelines, or if there are other considerations to be taken into account, please refer to the page for information on Custom Helmets or call for assistance.

HEADGEAR SELECTION

Some suggested considerations to use in helping to determine what type of helmet would be suitable for your client:

- Degree of physical activity.
- Potential for striking surfaces (floor, walls, furniture, etc.).
- Nature of self-abusiveness.
- Severity & frequency of seizures or spasticity.
- Length of time worn.

IMPORTANT:

While we believe that headgear can substantially reduce the risk of head injury, we emphasize that there is no known device for preventing all injuries. Only you the customer are in a position to evaluate the particularized needs of the individual client in each environment. It is essential that professional and medical expertise be utilized in rendering decisions regarding headgear use.

Measurements are in inches