



SwabTek Cocaine Test Kit

ZNTK-COC-25-00

ITEM #	ZNTK-COC-25-00
DESCRIPTION	<p>SwabTek's Cocaine Test Kit is designed to test for the presumptive presence of cocaine in any sample. Since the kit uses a treated cotton swab, it can be used to test for cocaine residue on any surface without the presence of a visible powder or rock. This allows the user to test any surface suspected to have come in contact with the drug, including tabletops, pipes, baggies and hands. The test is simple and intuitive — and contains no acids, carcinogens or toxic liquids — eliminating any potential risk to the user during testing. Samples may appear white or off-white in color and consist of a fine to coarse powder, a hard rock or opaque crystal, a sticky white paste, cocaine impregnated into paper sheets or dissolved as liquids and vaping juices.</p> <p>Screens for cocaine compounds and products, including:</p> <ul style="list-style-type: none"> • Cocaine, Cocaine Hydrochloride, Cocaine Freebase, Crack Cocaine, Speedball (Heroin/Cocaine mix)
SIZE	Case of 25 tests: 9.5" L x 5.75" W x 5" H
COMPONENTS	<p>SwabTek's Cocaine Test Kit consists of three pieces</p> <ul style="list-style-type: none"> • 1 x pre-treated swab • 1 x reagent-treated dry test strip • 1 x test instruction card
PACKAGING	25 packaged test kits/case
WEIGHT	Case of 25 tests: 1.12 lbs.
PRIMARY COUNTRY OF ORIGIN	United States
BRAND	SwabTek
CASE PACK	MASTER CASE: 20 bx/cs, 500ea
PALLET QUANTITY	24 cs/pallet
NOTES	<p>Veriteque USA Inc. (SwabTek) field tests are presumptive only and, as such, they indicate the presumed presence of chemical groups and precursors which may be present in a given sample. ALL SWABTEK TEST RESULTS SHOULD BE CONFIRMED BY AN APPROVED ANALYTICAL LABORATORY. All SwabTek tests must be administered in strict accordance with the specific instruction and reference materials that accompany the products for best results.</p>

View: [User Manual](#) | View: [Data Sheet](#)