



SwabTek Heroin Test Kit

ZNTK-HER-25-00

ITEM #	ZNTK-HER-25-00
DESCRIPTION	<p>SwabTek’s Heroin Test Kit is designed to test for the presumptive presence of heroin in any sample. Since the kit uses a treated cotton swab, it can be used to test for heroin residue on any surface without the presence of a visible powder. This allows the user to test any surface that is suspected to have come in contact with the drug, including spoons, needles, baggies and hands. The test is simple and intuitive — and contains no acids, carcinogens or toxic liquids — eliminating any potential risk to the user during testing. Samples may appear white, brown or black in color and consist of a fine to coarse powder, a hard rock or a sticky residue or in pastes, liquids, oils and vaping juices.</p> <p>Screens for heroine compounds and products, including:</p> <ul style="list-style-type: none"> • Heroin Free Base • Black Tar • Afghan Brown • Mexican Brown • China White • China Girl • Brown Sugar • Hell Dust • Speedball
SIZE	Case of 25 tests: 9.5" L x 5.75" W x 5" H
COMPONENTS	<p>SwabTek’s Heroine Test Kit consists of three pieces</p> <ul style="list-style-type: none"> • 1 x pre-treated swab • 1 x reagent-treated dry test strip • 1 x test instruction card
PACKAGING	25 packaged test kits/case
WEIGHT	Case of 25 tests: 1.12 lbs.
PRIMARY COUNTRY OF ORIGIN	United States
BRAND	SwabTek
CASE PACK	MASTER CASE: 20 bx/cs, 500ea
PALLET QUANTITY	24 cs/pallet
NOTES	<p>Veriteque USA Inc. (SwabTek) field tests are presumptive only and, as such, they indicate the presumed presence of chemical groups and precursors which may be present in a given sample. ALL SWABTEK TEST RESULTS SHOULD BE CONFIRMED BY AN APPROVED ANALYTICAL LABORATORY. All SwabTek tests must be administered in strict accordance with the specific instruction and reference materials that accompany the products for best results.</p>